
符合GMP要求的气流粉碎机组 Jet Mill unit certificated by GMP

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PRINCIPLE

The raw materials of traditional Chinese medicine are initially processed from the knife type coarse crusher, then followed by secondary impact and fine crushing, and finally fed into the jet mill for GMP micronized pulverization. The high-speed airflow impact arising from several nozzles in the pneumatic pulverizing chamber enables the materials to collide and abrade each other, fracture instantly and achieve the purpose of micronized pulverization. After being sorted from the classifier, the conformable products will be discharged with airflow and then fed into the cyclone separator and collected by the collector. The coarse particles are fed back to the coarse crusher for continuous pulverization and the gas will be exhausted from the induced draft fan.

FEATURES

- This unit is classified into two categories, normal temperature pulverization technology. The 0.5~10 micron powder particles (500~10,000nanometers) can be produced

as a result of extremely high airflow speed and tremendous impact force.

- Classifying devices are available inside the pulverizer, through which the coarse particles from the processing materials can be cyclically pulverized to produce finished products with uniform grain fineness and small range of particle diameters.
- This unit is so designed and structured under the appeal of GMP. The airflow is extremely pure with filtering process. Compact internal structure to perform the closed circuit milling. From raw materials to continuous production of finished products, the pulverization requires very short time but yields higher efficiency and continuous operation.


APPLICATON SCOPE

This unit applies to micronized pulverization of a wide variety of traditional Chinese medicines, western medicines, and health food.

符合GMP要求的气流粉碎机组主要参数 Major Technical Parameters of Jet Mill unit certificated by GMP

参数 Parameter	型号 Model	QYF-100	QYF-150	QYF-260	QYF-400	QYF-600
生产能力(kg/h) Capacity(kg/h)		0.5~10	5~100	50~200	80~380	200~500
空气耗量(m³/min) Air Consumption(m³/min)		2	3	6	10	20
工作压力(Mpa) Working Pressure(Mpa)		0.75~0.85	0.8~1.5	0.8~1.5	0.8~1.5	0.8~1.5
进料粒径(目) Feed Diameter(mesh)		100~325	60~120	60~120	60~120	60~120
粉碎细度(μm) Grinding Size(μm)		0.5~30	0.5~30	0.5~30	0.5~30	0.5~30
装机功率(kw) Energy Consumption Power(kw)		20	40	60	95	188

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